

## SMARTeZ™ Pump • eZS™ Pump Non-electrically Driven Portable Infusion Pumps



### Factors Affecting Flow Rate

The SMARTeZ™ Pump has a specially designed multi-layered balloon-like reservoir to be filled with the drug or fluid intended for infusion.

It exerts mechanical pressure thereby administering the contents through an orifice tube at a predetermined flow rate. The entire unit is sterile and is intended for single-use only.


Elastomeric pumps are calibrated to specific operating conditions.

When filled to the nominal volume, flow accuracy is within +/- 15% of the nominal (label) flow rate.

Several Factors may affect  
Flow Rate

- Fill Volume
- Temperature
- Pump Position
- Viscosity
- Storage Time
- External Pressure







1. Fill port
2. Outer soft cover
3. Multi-layered elastomeric membrane
4. ON-OFF clamp
5. PVC administration tubing
6. In-line air and particulate eliminating filter
7. Flow restrictor
8. Patient connector
9. Patient end cap
10. Fill port cap (tethered)
11. Labelling - Fill volume & infusion duration
12. Labelling - Flow rate



## SMARTeZ™ Pump • eZS™ Pump

### Non-electrically Driven Portable Infusion Pumps

When filled to nominal volume, flow accuracy is within +/- 15% of the nominal (label) flow rate when infusion is started within 8 hours after fill.

Factor	Impact on Flow Rate	Guidelines
Fill Volume		
	<p>Flow Rate accuracy specifications is based on filling the pump to the nominal volume.</p>	<p>Follow the Delivery Times Table in the Product User Information.</p> <ul style="list-style-type: none"> <li>- Always fill within the Minimum and Maximum limits specified for the pump.</li> </ul>
Temperature		
	<p>Temperature affects fluid viscosity, hence Flow Rate.</p> <p>Flow Rate increases 1.5% for every 1°F increase in temperature.</p> <p>Flow Rate decreases 1.5% for every 1°F decrease in temperature.</p>	<ul style="list-style-type: none"> <li>- If pump is refrigerated, allow it to come to room temperature before use.</li> <li>- The flow restrictor which is embedded in the patient connector should be kept close to the skin.</li> <li>- Do not place the pump under bedcovers or near a heat source as temperature may be too warm.</li> </ul>
Pump Position		
	<p>Pump position is best kept at same level as catheter site.</p> <p>Pump position affects flow rate - higher pump position results in faster flow rate while lower position decrease flow.</p>	<ul style="list-style-type: none"> <li>- During use it is best to put the pump in a carrier bag, pocket or bedside table.</li> <li>- Do not hang the pump on an IV pole.</li> </ul>
Viscosity		
	<p>The nominal flow rate is based on saline as the diluent.</p> <p>The use of viscous drugs or fluids as diluent will reduce flow rate.</p>	<ul style="list-style-type: none"> <li>- Use of 5% Dextrose will result in a 10% longer delivery time.</li> </ul>
Storage Time		
	<p>Flow Rate accuracy specification is based on starting infusion within 8 hours after filling.</p> <p>Your SMARTeZ Pump can be started anytime ONE hour after filling.</p>	<ul style="list-style-type: none"> <li>- Allow the pump to rest ONE hour after filling before starting infusion.</li> </ul>
External Pressure		
	<p>Flow Rate will increase if there are external pressures exerted on the pump.</p>	<ul style="list-style-type: none"> <li>- Do not squeeze or lie on the pump during infusion.</li> </ul>