

SMARTeZ Chemo Closed System Pumps

Type	Article No. (EU)	Volume (mL)	Flow Rate (mL/h)	Duration (Hrs)	Duration (Days)	Qty Per Carton	Color Code
Chemo CS with ProSeal	484460	60 mL	0.5 mL / h	120 Hrs	5 Days	10	Green
	484520	65 mL	0.5 mL / h	130 Hrs	5.5 Days	10	Green
	484500	270 mL	1.5 mL / h	180 Hrs	7.5 Days	10	Blue
	484450	50 mL	2 mL / h	25 Hrs	1 Day	10	Yellow
	484510	60 mL	2 mL / h	30 Hrs	1 Day	10	Yellow
	484400	100 mL	2 mL / h	50 Hrs	2 Days	10	Yellow
	484410	270 mL	2 mL / h	135 Hrs	5.5 Days	10	Yellow
	484420	120 mL	4 mL / h	30 Hrs	1 Day	10	Orange
	484530	240 mL	5 mL / h	48 Hrs	2 Days	10	Light Brown
	484430	270 mL	5 mL / h	54 Hrs	2 Days	10	Light Brown
	484470	240 mL	10 mL / h	24 Hrs	1 Day	10	Light Green
	484440	270 mL	10 mL / h	27 Hrs	1 Day	10	Light Green
	484480	120 mL	30 mL / h	4 Hrs	> 1 Day	10	Brown
	484490	100 mL	100 mL / h	1 Hr	> 1 Day	10	White