

SMARTeZ Pumps Rigid Shell

Type	Article No. (EU)	Volume (mL)	Flow Rate (mL/h)	Duration (Hrs)	Duration (Days)	Qty Per Carton	Color Code
Long Time Infusion - H	480150H	65 mL	0.5 mL / h	130 Hrs	5.5 Days	10	Green
	480680H	84 mL	0.5 mL / h	168 Hrs	7 Days	10	Green
	480300H	85 mL	0.5 mL / h	170 Hrs	7 Days	10	Green
	480180H	100 mL	0.5 mL / h	200 Hrs	8 Days	10	Green
	480700H	96 mL	1 mL / h	96 Hrs	4 Days	10	Pink
	480310H	100 mL	1 mL / h	100 Hrs	4 Days	10	Pink
	480380H	120 mL	1 mL / h	120 Hrs	5 Days	10	Pink
	480120H	100 mL	1.5 mL / h	67 Hrs	3 Days	10	Dark Blue
	480720H	250 mL	1.5 mL / h	167 Hrs	7 Days	10	Dark Blue
	480230H	270 mL	1.5 mL / h	180 Hrs	7.5 Days	10	Dark Blue
	480750H	300 mL	1.5 mL / h	200 Hrs	8 Days	10	Dark Blue
	480050H	60 mL	2 mL / h	30 Hrs	1 Days	10	Yellow
	480710H	96 mL	2 mL / h	48 Hrs	2 Days	10	Yellow
	480080H	100 mL	2 mL / h	50 Hrs	2 Days	10	Yellow
	480100H	120 mL	2 mL / h	60 Hrs	2 Days	10	Yellow
	480200H	240 mL	2 mL / h	120 Hrs	5 Days	10	Yellow
	480670H	240 mL	2 mL / h	120 Hrs	5 Days	10	Yellow
	480160H	270 mL	2 mL / h	135 Hrs	5 Days	10	Yellow
	480170H	300 mL	2 mL / h	150 Hrs	6 Days	10	Yellow
	480280H	60 mL	2.5 mL / h	24 Hrs	1 Days	10	Dark Brown
	480740H	300 mL	2.5 mL / h	120 Hrs	5 Days	10	Dark Brown
	480390H	150 mL	3 mL / h	50 Hrs	2 Days	10	Gold Yellow
	480690H	96 mL	4 mL / h	24 Hrs	1 Day	10	Orange
	480220H	100 mL	4 mL / h	25 Hrs	1 Day	10	Orange
	480060H	120 mL	4 mL / h	30 Hrs	1 Day	10	Orange
	480010H	60 mL	5 mL / h	12 Hrs	0.5 Days	10	Brown
	480020H	80 mL	5 mL / h	16 Hrs	0.5 Days	10	Brown
	480030H	125 mL	5 mL / h	25 Hrs	1 Day	10	Brown
	480400H	240 mL	5 mL / h	48 Hrs	2 Days	10	Brown
	480090H	270 mL	5 mL / h	54 Hrs	2 Days	10	Brown
	480350H	300 mL	5 mL / h	60 Hrs	3 Days	10	Brown
	480210H	100 mL	8 mL / h	12.5 Hrs	0.5 Days	10	White
480370H	250 mL	8 mL / h	31 Hrs	1 Day	10	Dark Gray	
480040H	270 mL	10 mL / h	27 Hrs	1 Day	10	Lime	
480360H	300 mL	10 mL / h	30 Days	1 Day	10	Maroon	
480730H	275 mL	11.5 mL / h	24 Hrs	1 Day	10	Peach	
Type	Article No. (EU)	Volume (mL)	Flow Rate (mL/h)	Duration (Hrs)	Duration (Mins)	Qty Per Carton	Color Code
Short Time Infusion - H	481190H	50 mL	100 mL / h	0.5 Hrs	30 Mins	10	White
	481010H	100 mL	200 mL / h	0.5 Hrs	30 Mins	10	Light Blue
	481030H	50 mL	50 mL / h	1 Hr	60 Mins	10	Dark Green
	481090H	100 mL	50 mL / h	2 Hr	120 Mins	10	Dark Green
	481040H	100 mL	100 mL / h	1 Hr	60 Mins	10	White
	481050H	250 mL	250 mL / h	1 Hr	60 Mins	10	Dark Gray

Type	Article No. (EU)	Volume (mL)	Flow Rate (mL/h)	Duration (Hrs)	Duration (Days)	Qty Per Carton	Color Code
Chemo Infusion - H	484070H	60 mL	0.5 mL / h	120 Hrs	5 Days	10	Green
	484130H	65 mL	0.5 mL / h	130 Hrs	5.5 Days	10	Green
	484110H	270 mL	1.5 mL / h	180 Hrs	7.5 Days	10	Dark Blue
	484060H	50 mL	2 mL / h	25 Hrs	1 Day	10	Yellow
	484120H	60 mL	2 mL / h	30 Hrs	1 Day	10	Yellow
	484010H	100 mL	2 mL / h	50 Hrs	2 Days	10	Yellow
	484020H	270 mL	2 mL / h	135 Hrs	5.5 Days	10	Yellow
	484200H	275 mL	2.2 mL / h	125 Hrs	5 Days	10	Yellow
	484220H	100 mL	4 mL / h	25 Hrs	1 Day	10	Orange
	484030H	120 mL	4 mL / h	30 Hrs	1 Day	10	Orange
	484210H	200 mL	4 mL / h	50 Hrs	2 Days	10	Orange
	484180H	125 mL	5 mL / h	25 Hrs	1 Day	10	Light Brown
	484040H	270 mL	5 mL / h	54 Hrs	2 Days	10	Light Brown
	484190H	275 mL	5.7 mL / h	48 Hrs	2 Days	10	Light Brown
	484080H	240 mL	10 mL / h	24 Hrs	1 Day	10	Light Green
	484050H	270 mL	10 mL / h	27 Hrs	1 Day	10	Light Green
	484100H	100 mL	100 mL / h	1 Hr	> 1 Day	10	White