

SMARTeZ Pumps - Bolus Rigid-Shell

Type	Article No. (EU)	Volume (mL)	Flow Rate (mL/h)	Duration (Hrs)	Duration (Days)	Qty Per Carton	Color Code
Bolus-H	386010H	100 mL	2 mL / h	50 Hrs	2 Days	10	Yellow
	386040H	100 mL	2 mL / h	50 Hrs	2 Days	10	Yellow
	386020H	270 mL	2 mL / h	135 Hrs	5.5 Days	10	Yellow
	386050H	270 mL	2 mL / h	135 Hrs	5.5 Days	10	Yellow