

# SMARTeZ Pumps - PCA Bolus & Flow

Type	Article No. (EU)	Volume (mL)	Bolus Volume (mL)	Flow Rate (mL/h)	Refill Time	Qty Per Carton
PCA-Bolus & Flow	485180	100	0.5 mL	0 - 1 - 2 - 3	15 Mins	10
	485190	120	1 mL	0 - 1 - 2 - 3	15 Mins	10
	485200	200	1 mL	0 - 2 - 4 - 6	15 Mins	10
	485220	xx	xx	0 - 2 - 4 - 6	xx	10
	485210	270	2 mL	0 - 3 - 6 - 9	15 Mins	10
	485010	100	2 mL	1 - 2 - 3 - 4	30 Mins	10
	485110	100	5 mL	1 - 2 - 3 - 4	30 Mins	10
	485120	270	5 mL	1 - 2 - 3 - 4	30 Mins	10
	485030	400	2 mL	1 - 2 - 3 - 4	30 Mins	10
	485060	400	2 mL	1 - 2 - 3 - 4	30 Mins	10
	485130	400	5 mL	1 - 2 - 3 - 4	30 Mins	10
	485140	100	5 mL	2 - 4 - 6 - 8	30 Mins	10
	485040	100	2 mL	2 - 4 - 6 - 8	60 Mins	10
	485090	270	2 mL	2 - 4 - 6 - 8	15 Mins	10
	485050	270	2 mL	2 - 4 - 6 - 8	30 Mins	10
	485150	270	5 mL	2 - 4 - 6 - 8	30 Mins	10
	485070	270	2 mL	2 - 4 - 6 - 8	60 Mins	10
	485100	400	2 mL	2 - 4 - 6 - 8	15 Mins	10
	485080	400	2 mL	2 - 4 - 6 - 8	30 Mins	10
	485160	400	5 mL	2 - 4 - 6 - 8	30 Mins	10
	485060	400	2 mL	2 - 4 - 6 - 8	60 Mins	10
	485170	400	2 mL	4 - 6 - 8 - 10	60 Mins	10