

SMARTeZ Pumps PainBloc

Type	Article No. (EU)	Volume (mL)	Flow Rate (mL/h)	Duration (Hrs)	Duration (Days)	Qty Per Carton	Color Code
PainBloc	483010	100 mL	1 mL / h	100 Hrs	4 Days	10	Pink
	483020	100 mL	2 mL / h	50 Hrs	2 Days	10	Yellow
	483030	270 mL	2 mL / h	135 Hrs	5.5 Days	10	Yellow
	483040	270 mL	4 mL / h	68 Hrs	3 Days	10	Orange
	483060	400 mL	4 mL / h	100 Hrs	4 Days	10	Orange
	483050	270 mL	5 mL / h	54 Hrs	2 Days	10	Brown
	483070	400 mL	5 mL / h	80 Hrs	3.5 Days	10	Brown
	483080	400 mL	6 mL / h	67 Hrs	3 Days	10	Old Rose
	483090	400 mL	8 mL / h	50 Hrs	2 Days	10	Purple
	483100	400 mL	10 mL / h	40 Hrs	2 Days	10	Light Green