

SMARTeZ Pumps - Bolus

Type	Article No. (EU)	Volume (mL)	Flow Rate (mL/h)	Duration (Hrs)	Duration (Days)	Refill Time	Qty Per Carton	Color Code
Bolus	486080	60	0			7 Mins	10	Dark Gray
	486070	60	0			15 Mins	10	Dark Gray
	486090	100	0			7 Mins	10	Dark Gray
	486100	100	0			15 Mins	10	Dark Gray
	486120	100	1	100 Hrs	4 Days	15 Mins	10	Pink
	486110	100	2	50 Hrs	2 Days	12 Mins	10	Yellow
	486160	100	2	50 Hrs	2 Days	15 Mins	10	Yellow
	486010	100	2	50 Hrs	2 Days	30 Mins	10	Yellow
	486040	100	2	50 Hrs	2 Days	60 Mins	10	Yellow
	486130	120	2	60 Hrs	2.5 Days	15 Mins	10	Yellow
	486140	200	2	100 Hrs	4 Days	15 Mins	10	Yellow
	486020	270	2	135 Hrs	5.5 Days	30 Mins	10	Yellow
	486050	270	2	135 Hrs	5.5 Days	60 Mins	10	Yellow
	486030	400	2	200 Hrs	8 Days	30 Mins	10	Yellow
	486060	400	2	200 Hrs	8 Days	60 Mins	10	Yellow
	486150	270	4	67 Hrs	2.5 Days	15 Mins	10	Orange